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EN USER MANUAL

ORIGINAL NOVAFON INTRAORAL ATTACHMENTS

- PREFACE**
- PRODUCT DESCRIPTION**
- APPLICATIONS & USE**
- IMPORTANT INSTRUCTIONS**
- CONTRAINDICATIONS**
- SIDE EFFECTS**
- CLEANING AND DISINFECTION**
- CONTACT**

1. PREFACE

The intraoral attachments (tongue depressor head, arrow head, spoon head, and pellet head) are used for intraoral application of vibrotactile stimulation in conjunction with NOVAFON sound wave devices. They are compatible with any NOVAFON device. Please read the following pages carefully to ensure proper use.

2. PRODUCT DESCRIPTION

2.1. COMPONENTS

2.1.1. Tongue depressor head [item no. 5003]

The tongue depressor head is modelled after a tongue depressor and has two different surface structures. The special grooves and ripples ensure that the tongue cannot slip during application.

2.1.2. Arrow head [item no. 5004]

The arrow head has three different surface structures (rippling, nubs, smooth) and its shape is ideally suited for intra- and extraoral applications.

2.1.3. Spoon head [item no. 5005]

The spoon head has a spoon-like trough with nubs on the bottom side. It can be used both for vibrotactile stimulation and for the application of small amounts of food.

2.1.4. Pellet head [item no. 5006]

The pointed pellet head with a small ball at the end of the attachment can be used for very targeted and precise intra- and extraoral stimulation.

2.1.5. Adapter [item no. 5007]

It is already mounted to the attachment upon delivery. It is available as a spare part and fits on all NOVAFON sound wave devices as well as on any intraoral attachment.

2.2. MATERIAL INFORMATION

The intraoral attachments are made of thermoplastic elastomers, a high-quality plastic specifically designed for medical applications. The material has no added of phthalates or latex. Please store it in a cool and dry place.

2.3. SHELF LIFE

The shelf life of the intraoral attachments will vary depending on the frequency of use. We recommend replacing the attachments after one year of use. Check the attachments before each application. Replace damaged, deformed or heavily discoloured products.

3. APPLICATION

First connect the intraoral attachment to the adapter. Then screw the adapter onto your NOVAFON sound wave device. If you wish to change or remove the attachment, simply unscrew the adapter with the attachment from the device. Then remove the attachment from the adapter by carefully pulling it off. Please note that it may be necessary to apply light counter-pressure to the attachment so that it does not detach from the fixture. Attention! Before using the intraoral attachments, please pay attention to the contraindications listed in chapter 5.

Application examples for intraoral attachments:

Tongue depressor head	Arrow head	Spoon head	Pellet head
Flat tapping on the tongue	Ability to sense different surfaces	Training of tongue tip motility	Indicates place of articulation
Tongue edge and cheek stimulation	Stimulation of tongue, cheeks and lips	Stimulation of the back of the tongue	Stimulation of the palatine arches
Holding the tongue depressor between the lips	Tongue furrow formation	Application of small amounts of food	Central tongue furrow formation

Caution! Be careful to avoid touching the teeth with the attachments as much as possible, as this can be particularly unpleasant for existing tooth infections.

4. IMPORTANT INSTRUCTIONS

4.1. HANDLING OF INTRAORAL ATTACHMENTS AND NOVAFON SOUND WAVE DEVICES

- Always connect the intraoral attachment to the adapter first and then screw it onto your NOVAFON sound wave device.
- If you wish to change or remove the attachment, always unscrew the adapter with the attachment from the device before you detach the attachment from the adapter. Otherwise, the device may be damaged.
- The attachments can be rotated within the adapter. If you want to rotate your attachment inside the adapter while it is already on your NOVAFON sound wave device, hold the adapter firmly with one hand while you rotate the attachment. Otherwise, the device may be damaged.
- When you exert pressure while using the attachment, always press with one finger against the attachment, so that no pressure is exerted on the head of the device. Otherwise, the device may be damaged.
- Before using the attachments, please also observe all information from the operating instructions of your NOVAFON sound wave device.
- Attention! Improper use may result in injury.

4.2. INTENDED USERS

The intraoral attachments are intended for professional use only. They are only to be used by trained, medical personnel or under strict supervision of such a person.

4.3. HYGIENE

Before using the intraoral attachments for the first time, they must undergo a specific hygienic treatment (according to chapter 7.) For use on several people, the attachments must be prepared hygienically after each application (according to chapter 7.). Please wear gloves for each application.

5. CONTRAINDICATIONS

An application of the intraoral attachments is not permissible for:

- Inflammation and injury to the oral mucosa (aphthae (canker sores), thrush, bleeding, mucositis, etc.)
- Open wounds in the oral cavity
- Tumours in the oral cavity
- Bite reactions / excessive reflexes

6. SIDE EFFECTS

Intraoral stimulation may occasionally result in the following side effects:

- Triggering the bite reflex
- Increased muscle tone
- Increased salivation
- Defensive behaviour

7. CLEANING AND DISINFECTION

PROCEDURE	
<ul style="list-style-type: none">Cleaning and disinfection: manual Cleaning and disinfection: by machine	
Warning information	The attachments must be cleaned and, if necessary, disinfected before the first application. For use on several people, the products must be cleaned and disinfected after each application.
Restriction regarding reprocessing	No special requirements.
INSTRUCTIONS	
Application location	No special requirements.
Storage and transport	It is highly recommended to transport contaminated intraoral attachments in a closed container. It is recommended to carry out the reprocessing of intraoral attachments as soon as possible, within a maximum of 2 hours after use.
Preparation	Wear personal protective equipment (gloves, water repellent lab coat, face mask, or safety glasses). Disassemble all detachable intraoral attachments into their individual parts.
Pretreatment	Equipment: Plastic brush (e.g., Interlock, #09098 and Interlock, #09050), running water (at least drinking water quality) Procedure: Brush coarse surface contamination on the attachments and the adapter (thread and cavity) under running water.

MANUAL

Manual cleaning	Equipment: Cleaning agent [multistage enzymatic cleaner, e.g. ECOLAB Sekusept® MultiEnzyme P [ecolab.com]], plastic brush (e.g. Interlock, #09098 and Interlock, #09050), running water (at least drinking water quality), tray for cleaning agent. <ol style="list-style-type: none">Apply the cleaning solution according to the manufacturer’s instructions (e.g. 2% Sekusept® MultiEnzyme P). Submerge all parts of the attachments completely. Remove surface contamination with the plastic brush, by scrubbing all surfaces of the attachment and the adapter. Treatment time in the cleaning solution according to the manufacturer’s instructions of the cleaning agent (e.g. 2% Sekusept® MultiEnzyme P, at least 15 min). Remove all parts from the cleaning solution. Rinse each part of the attachments under running water for 30 sec. Check for cleanliness and, if there is still visible contamination, repeat the pretreatment and the aforementioned manual steps.
-----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Manual disinfection	Equipment: Virucidal instrument disinfectant (VAH or IHO listed) based on peroxide compounds or quaternary compounds and alkyl-amine/alkylamine derivatives (e.g. ECOLAB Sekusept® aktiv or Dr. Schumacher, Perfektan active), disinfectant tray <ol style="list-style-type: none">Apply the disinfectant solution according to the disinfectant manufacturer’s instructions (e.g. Sekusept® aktiv, 2% / 2% Perfektan active). Submerge all parts of the attachments completely. Treatment time in the disinfectant solution according to the manufacturer’s instructions for the disinfectant (e.g. at least 5 minutes for Sekusept® aktiv and at least 30 minutes for Perfektan active). Remove all parts of the attachment from the disinfectant solution. Rinse each item for 1 minute under running water. Wipe all parts with a lint-free disposable cloth for drying and, if necessary, dry it with medical compressed air.
---------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

BY MACHINE				
Cleaning and disinfection: by machine	Note: Coarse surface contamination on the attachments must be removed before preparation by machine (see pretreatment: manual). Cleaning and disinfecting device according to DIN EN ISO 15883-1+2 with thermal program (temperature 90 °C to 95 °C). - Mildly alkaline cleaner (e.g. ECOLAB Sekumatic® MultiClean) <ol style="list-style-type: none">Place all the individual parts of the attachments in a suitable strainer bowl or on the load carrier in such a way that all inner and outer surfaces are cleaned and disinfected. Close the RDG and start the programme; for programme sequence see table below. At the end of the programme, remove the products. Check the cleanliness and dryness of the load (all parts) and, if necessary, dry with medical compressed air. For residual contamination, carry out the manual pre-treatment and manual cleaning.			

Programme sequence	Water	Dosage	Time	Temperature
Prewash	cold water		5 min.	
Cleaner dosage		According to manufacturer		According to manufacturer
Cleaning	completely desalinated water		10 min.	55 °C
Wash	completely desalinated water		2 min.	
Disinfecting	completely desalinated water			A0 value > 3,000' [e.g. 90 °C, 5 min.]
Drying			15 min.	Up to 120 °C

¹ Authorities may adopt other implementing provisions (parameters for disinfection performance) within their jurisdiction.

Maintenance, inspection, and testing	Check all parts of the attachments for damage and wear. Damaged medical devices may no longer be used and must be sorted out.
--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------

Storage	Storage and storage period according to the specifications of the user.
---------	-------------------------------------------------------------------------

8. CONTACT

If you have any questions, please contact your local dealer or the manufacturer.

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