novafon®

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English

1. Key for overview	27
2. Preface	27
3. Product description	27
3.1. Scope of delivery	27
3.2. Device type, technical specifications	27
3.3. Accessories	29
3.3.1. Disc attachment	29
3.3.2. Ball attachment	29
3.3.3. Charging stand	29
3.3.4. Optional accessories	29
3.3.5. Combination with other products	29
4. Operation	30
4.1. The Display	30
4.2. Initial start-up and frequency	30
4.3. Setting the intensity / strength	30
4.4. Charging	31
4.5. Attachments	31
5. Application & use & user groups	32
5.1. Intended use	32
5.2. General information for use	32
6. Contraindications and Side Effects	34
7. Technical specifications	35
7.1. Symbols used	35
7.2. Technical specifications	35
8. Important Safety Instructions	37
8.1. Warnings	37
8.2. Other precautions	41
9. Maintenance and repair	43
9.1. Cleaning and disinfection	43
9.2. Troubleshooting	43
9.3. Safety check, maintenance and repair	44
9.4. Disposal	44
10. Warranty	45
11. Customer service	45
12. Simplified EU Declaration of Conformity	46

1. Key for overview

On page 1 you will find an overview of your device with the following important elements:

1. Display Intensity
2. Multi-button
top: intensity setting in three levels
bottom: frequency setting
Standby / ON / 100 Hz / 50 Hz / OFF
3. Frequency indicator
1 illuminated dot (●): 100 Hz
2 illuminated dots (●●): 50 Hz
4. Battery / charging status indicator
5. USB charging port
6. Correct positioning of the attachment on the device

2. Preface

Congratulations on your purchase of a novafon vibration sound wave device. Our devices have delivered excellent performance and reliability for decades. Please read the following instructions carefully to ensure that you can get the most out of your novafon vibration sound wave device and enjoy using it for a long time.

Please keep these instructions for future reference.

3. Product description

3.1. Scope of delivery

Delivery includes a battery-powered novafon vibration sound wave device as well as a disc attachment, a ball attachment, a power plug, a charging cable and instructions for use.

3.2. Device type, technical specifications

novafon sound wave devices have been used for decades to relieve pain and other ailments. Vibration massage can have a relaxing or invigorating effect and contribute to a sense of well-being. The transmission of gentle vibrations to the skin's surface,

muscles, cartilage and tissue produces a stimulating effect that is often still noticeable for several minutes or even hours after treatment.

Thanks to the two treatment frequencies and three intensity settings, treatment can be individually adapted to the body part and underlying symptoms in question.

The following table describes the characteristics of your new novafon device

Device type	novafon one
Multifunction buttons	ON/OFF (⏻) and (⏷)
frequency settings (⏻)	50 Hz, 100 Hz, Standby, Off
Intensity / strength setting (⏻)	low (⏻), medium (⏻), high (⏻)
Intensity	●●○○○
Depth of penetration	6 cm +
Colour	white
Lithium ion battery	Capacity 2400 mAh, voltage 3,7 V, 8,8 Wh, Type INR18650 with protection board
Charging time	approx. 3 h
Battery life	up to 3.5 h depending on frequency and intensity

For more technical specifications please see the technical data sheet in the download section of our website.

3.3. Accessories

novafon vibration sound wave devices are equipped with two interchangeable attachments as standard. NOVAFON GmbH also offers other attachments that are not included in this set. Please contact your dealer or visit www.novafon.com for more information. novafon devices may only be used with original novafon accessories.

3.3.1. Disc attachment

The disc attachment is used in the treatment of larger body parts that require broader strokes. It is supplied as standard and also available as a spare part.

3.3.2. Ball attachment

The ball attachment can be used for smaller areas and targeted treatment. It is supplied as standard and also available as a spare part.

3.3.3. Charging stand

The charging stand can be purchased together with the device or separately. The charging stand offers convenient storage for your device during and after charging. These instructions for use include all the information you need on how to use and clean your stand, as well as any necessary precautionary measures.

3.3.4. Optional accessories

On our website you will find additional attachments that simplify treatments in different areas of the body.

3.3.5. Combination with other products

novafon devices are not intended to be used with products other than novafon original accessories. **The use of accessories from other manufacturers can impair the function and performance of novafon devices.**

4. Operation

4.1. The Display

The current operating status is displayed by the 3 arcs above the button and the 2 points below the button. If none of the indicators are lit, the device is off/stationary or the battery is completely empty. The two dots show the currently selected frequency (1 point = 100 Hz, 2 dots = 50Hz). The three arcs the currently selected intensity (low, medium, high). The single dot on the bottom end of the handle flashes when the battery level is low or when the device is charging.

4.2. Initial start-up and frequency

Fully charge the device prior to switching it on for the first time (chapter 4.4).

The device starts at 100 Hz (♫) by default. Pressing the button again decreases the frequency to 50 Hz (♫). The current frequency selection is indicated by the LED display (frequencies). Press the button a third time to switch off the device.

4.3. Setting the intensity / strength

You can regulate the intensity of your device by repeatedly pressing the intensity button. Here you can choose from three intensity settings: low (☹), medium (☺), and high (☺). You can recognize the intensity selection by the LED display (intensity).

The device is set to start on the lowest intensity setting by default (☹). This also corresponds to novafon's recommendations for treatment i.e. start at low intensity and increase as required, according to each individual's sensitivity.

Note: The best vibration is achieved when the device is operated without additional contact pressure. If your device automatically switches to a lower setting, or you are unable to increase the intensity, this is probably due to a low battery level. In this case, you should charge your device.

4.4. Charging

To charge the device, you should only use the supplied power plug and cable, or the novafon charging stand if included (depending on the model, the charging stand may not be included in delivery). To charge your device, first connect the cable to the device or charging stand and then to the power plug suitable for your power connection. Only then should you plug the power plug into the socket. To charge the device without a charging stand, use the charging cable included in the device packaging. Charge the device until the battery is fully charged. Do not use the novafon while it is charging. The device has a low battery indicator.

When using the charging stand to charge your device, make sure you are careful when placing the device onto the connector and removing it. To remove, hold the charging stand and gently release the device. The device is charging when the battery indicator appears (see view 1).

Avoid twisting the device in the charging station in order to prevent damage.

To avoid malfunctions, the device must be completely dry and free from residues of cleaning agents, disinfectants or oils / creams.

Note: When the device is fully charged, disconnect the power plug and device from the electrical supply to avoid damage. If the device is connected to the charging accessory, all pole breaking and complete disconnection from the mains can only take place if the mains plug is removed from the socket.

4.5. Attachments

The attachments are fixed to the device via a plug connection. They can easily be detached by turning them a quarter turn and pulling them gently off the

device. The attachments click simply onto the device. The attachments are in the correct position when the brackets on the device and the attachment lie seamlessly on top of each other. The magnets help you to find the final position. Before starting treatment, make sure that the respective attachment is engaged in the intended final position. Please refer to the drawing on page 1 for the correct position of the attachment on the device.

5. Application & use & user groups

5.1. Intended use

The novafon sound wave device is designed to create a feeling of well-being. The gentle stimulation of the vibration massage and the associated heat effect create a physiological effect that, for many users, relaxes the muscles in the contact area and can also temporarily relieve minor pain. This creates a feeling of well-being.

Note: When using the novafon one on sensitive parts of the body, such as the neck, face and hands, make sure to only work with low intensities.

5.2. General information for use

Treatment should not exceed a maximum of 20 minutes per session, nor a maximum of 3 sessions per day. Please note that after 20 minutes of use, the novafon should be switched off for 15 minutes in order to avoid overheating. If used several times a day on one person, allow 1-2 rest days per week without application. When treating severely painful or sensitive parts of the body (e.g. face, neck, head), you should reduce / adjust the duration of treatment according to each individual's sensitivity.

Note: Application of liniments: Apply only a thin layer of creams, ointments or cosmetic / medicinal liniments to the skin, or place a fine cloth between the device and the skin. Do not allow liquids or residues

to penetrate the novafon as this may cause damage. After using agents of this type, make sure you clean the device thoroughly.

Note: Discontinue treatment immediately if side effects occur (chapter 6).

The novafon appliance can be used directly on the skin or through fabric / clothing. Slowly pass the novafon over the affected area in a circular or sweeping motion. You can also use the device selectively and let it rest in place (for up to 60 seconds or until it becomes uncomfortable) to achieve a feeling of relief from pain and discomfort. Please note that you should only pass the device over the affected area with its own weight and not apply any additional pressure, as this will reduce the vibration effect. When applying to sensitive areas (e.g. face, neck, head) or multiple body parts in one session, the application time should be shortened and adjusted to individual needs. The novafon can be used on all parts of the body. Avoid contact with the temple areas or the eyes. **Do not use in the genital area.**

The session should feel comfortable and relaxing at all times. Stop using or change the position of the novafon if you find it painful or uncomfortable. **Do not use the novafon if any of the contraindications apply to you.**

It is normal for the attachments and even the device head to become warm during use. This heat is caused by the mechanics in the device. If the device head becomes hot, switch off the device and allow it to cool down for at least 15 minutes before using it again. If the device remains hot, stop using it immediately and either discard it or contact your dealer.

6. Contraindications and Side Effects

Treatments with the novafon vibration sound wave devices are not permitted for patients with:

- Open wounds/eczema (in the body region to be treated)
- Arteriosclerosis (in the body region being treated)
- Current episodes of epilepsy
- Pregnancy
- Implants (in the body region being treated)
- Heart pacemakers, defibrillators and similar
- Thromboses (in the body region being treated)
- Brain pacemakers and deep brain simulators
- Cardiac arrhythmia
- Tumours, cancerous and not cancerous
- Acute inflammations (in the body region being treated)
- Acute episodes of inflammatory diseases

Side effects:

Please note that the following side effects may occur:

- Increased pain/unpleasant sensation (occasionally)
- Discomfort/dizziness/restlessness (occasionally)
- Skin irritation (rarely)
- Bruises (very rarely)
- Reflexes/involuntary muscle contractions (rarely)
- Muscle soreness (rarely)

Note: Pain treatment in particular can lead to an initial worsening of the symptoms before you notice any improvement.

CAUTION: Discontinue treatment immediately if you experience any side effects. Consult a physician if symptoms persist.

Note: Slight skin redness and warmth are not a cause for concern. This is a desired effect of local vibration therapy, caused by increased blood flow.

Do not use to treat pain or swelling of unknown origin. In this case, consult a physician for evaluation and diagnosis.

7. Technical specifications

7.1. Symbols used

- ~ Alternating current (AC)
- = Direct current (DC)
-  Device with protection class II according to DIN EN 61140
- Hz Hertz
- KB Short-term operation
- mA Milliampères
- mAh Milliampere-hour
- A Amp
- V Volt
-  Disposal according to the EU Directive on Waste Electrical and Electronic Equipment (WEEE)
-  Displays the manufacturer of the device.
-  Displays the serial number as defined by the manufacturer to identify a specific device.
-  Protection class of the application part: type BF
-  Temperature limit: Refers to the temperature limits to which the device can be safely exposed
-  Humidity limit
-  Atmospheric pressure limit
-  Refer to the instructions
-  Use inside only
-  General warning symbols
-  Intensity button (3 intensities)
-  Frequency setting
- Stand-by / ON / 100 Hz / 50 Hz / OFF

7.2. Technical specifications

- Operating conditions: 5 °C (41 °F) to 30 °C (86 °F);
15 % to 93 % relative humidity, non-condensing;
700 hPa to 1060 hPa air pressure
- Storage / transport: -5 °C (23 °F) to 40 °C (104 °F);
15 % to 93 % relative humidity, non-condensing;
500 hPa to 1060 hPa air pressure

- Operating mode: Short-term operation 20/15 min
- Operating noise: maximum 70 dBA
Operating temperature: maximum 50 °C (122 °F) for maximum service life.
- Power plug Model 1180004-040
- Input power: AC 100–240 V, 50–60Hz, 300 mA
- Output power: DC 5 V, 2 A
- LOT: Shows the date of manufacture of the power supply
- A type F (EU plug) power plug is included in the standard scope of delivery. Additional mains adapters with different power plugs are optionally available or already included in the corresponding country sets.
- USB-A charging cable to USB-C, cable length: approx. 1 m
- One charge of a new battery guarantees at least six 20-minute treatment sessions at medium intensity. Note that frequent use at the highest intensity level can decrease battery life

Battery: Lithium-ion battery type INR18650, permanently installed, lifespan: 800 charging cycles

- Battery specifications: Capacity 2400 mAh, Voltage 3.7 V, 8.8 Wh; permanently installed, lifespan: 800 charging cycles
- Protection class IP22: protected against finger contact; protected against dripping water falling at an angle (up to 15°)

Due to the IP rating and safety requirements for medical devices, it is not possible for the battery to be replaced by the end user.

- Devices comply with IEC 60601-1-2:2014/AMD1:2020 and/or CISPR 11 Emission Class B and Group 1. Devices have been tested for home use in accordance with interference immunity testing levels.

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8. Important Safety Instructions

Before starting to use the device, please read the user manual and these safety instructions carefully and strictly follow them.

8.1. Warnings

When using an electrical device, basic precautions should always be followed including the following to reduce the risk of burns, fires, electric shock, possibly severe injury and possible death to the user or other people:

- ⚠ WARNING: Do not use the novafon device during charging.
- ⚠ WARNING: Always store the device in a dry and clean location.
- ⚠ WARNING: Before connecting the power plug to the mains power source, check that the mains voltage corresponds to the voltage specification on the power plug. Connect the novafon device only to a power outlet that delivers the correct voltage.
- ⚠ WARNING: To charge the device, use only the charging cable and power supply included in the delivery or the separately available charging stand. Use the correct connector plug, if applicable. Non-observance can lead to increased electromagnetic emissions or reduced electromagnetic immunity and thus to malfunctions or unexpected operational behaviour.
- ⚠ WARNING: Always unplug the device and all charging accessories from the electrical outlet after charging and before cleaning. Make sure that the mains plug is easily accessible at all times so that you can easily disconnect the device from the mains in an emergency. Do not place heavy objects or furniture in front of the

- power adapter, even when using the charging stand.
- ⚠ WARNING: Ensure that devices are completely dry and with no cleaning agent, disinfectant or oil / cream residues before using the charging stand or charging cable.
 - ⚠ WARNING: Do not bring the device and charging accessories (plug, charging cable, charging stand) into contact with hot surfaces (e.g. heater, hotplate, etc.). The housing of the device or adapter or the insulation of the cable could be damaged.
 - ⚠ WARNING: Do not use the device when taking a shower or bath.
 - ⚠ WARNING: Do not use the device where aerosol (spray) products are used, in other explosive atmospheres, or anywhere where oxygen is administered.
 - ⚠ WARNING: Do not place or drop the device into water or other liquids.
 - ⚠ WARNING: Do not reach for a device that has fallen into water, especially if the charging cable is plugged into a socket. Switch off the power immediately by switching off the fuse for the relevant circuit.
 - ⚠ WARNING: Do not place or store the device where it could fall or be pulled into a tub, pool or sink.
 - ⚠ WARNING: Protect the device, charger and charging stand from drop impact and shocks.
 - ⚠ WARNING: This device may not be modified without the manufacturer's permission.
 - ⚠ WARNING: novafon devices may only be used with original novafon accessories. The use of accessories from other manufacturers can impair the function and performance of novafon devices. This also applies to the power adapter and charging cables.
 - ⚠ WARNING: Do not use the device and either discard or return to your dealer for examination if any of the following conditions occur:
 - if the device fails to work correctly
 - if the device is damaged, e.g. after having been dropped
 - if the device was immersed in water or other liquids
 - if the device exhibits other safety-related problems, such as emitting smoke, unusual or excessive heat etc.
 - ⚠ WARNING: In order to prevent adverse events to the patient and operator due to electromagnetic disturbances, the novafon shall not be operated outside its intended EMC environment. Furthermore, the novafon shall not be operated if the enclosure, cables or measures for electromagnetic shielding are damaged.
 - ⚠ WARNING: Do not use the power plug, charging cable or charging stand:
 - in the event of visible damage or after dropping the power plug, charging cable or charging stand
 - if the power plug, charging cable or charging stand has been exposed to extreme moisture
 - if liquid has penetrated the unit
 - ⚠ WARNING: Never charge the device for more than 24 hours at a time. Do not leave the device unattended. Disconnect the device from the power supply at night. This also applies when using the charging stand.
 - ⚠ WARNING: Do not place the device on flammable surfaces when charging.
 - ⚠ WARNING: Never leave the device unattended when it is turned on.
 - ⚠ WARNING: Do not clamp or pinch any part of your body, especially your fingers, between the attachment and the device. During the treatment, hold the device only by the device handle.

- ⚠ WARNING: Never drop or insert any foreign object into any opening on this device. Do not insert any objects (including wires) into openings in the device.
- ⚠ WARNING: Wireless communication devices (such as Wi-Fi stations, mobile phones, cordless phones, and related base stations) can affect the device and should be kept at least 0.5 m (2 feet) away from the device during use. This also applies to cables and plug-in power supplies. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- ⚠ WARNING: Careful supervision is necessary when used on or by children, immobile or disabled individuals. novafoam devices should only be used by users from the age of 13 years who have read the instructions and are competent to use the device safely.
- ⚠ WARNING: Do not use on sleeping or unconscious individuals.
- ⚠ WARNING: Keep the device out of reach of children when not in use.
- ⚠ WARNING: Small parts (e.g. cover caps) could come loose and be swallowed. Risk of suffocation!
- ⚠ WARNING: Keep children away from packaging materials. There is a risk of suffocation.
- ⚠ WARNING: Do not use on genitals, this device is not intended as an intimate pleasure toy.
- ⚠ WARNING: Do not use outdoors or in wet or damp areas.
- ⚠ WARNING: Never operate the power plug or charging stand in a high-humidity environment or in the presence of moisture. Please note that

moisture is also generated by e.g. nebulisers and paraffin baths.

- ⚠ WARNING: Do not use or charge the device under blankets, pillows, towels, clothing, etc. or in places where ventilation is restricted, as it may overheat and could cause fire, electric shock or other personal injury.
- ⚠ WARNING: Use the device only for its intended use as described in this manual.
- ⚠ WARNING: Avoid contact with eyes.
- ⚠ WARNING: Avoid contact with mucous membranes unless using specially designed attachments.
- ⚠ WARNING: Do not use the device together with other devices and do not stack the devices on such products, as this may affect the performance and function of the device as well as the other product.
- ⚠ WARNING: Avoid overheating the device. Do not use the device for more than 20 minutes per application. Allow the device to cool down for at least 15 minutes before switching it on again.
- ⚠ WARNING: For technical reasons, the device can reach a temperature of up to 50 °C, especially on the head. Therefore, do not hold the device by your head for long periods of time. Change your grip position if the device feels uncomfortably warm. Switch off the device if any part of it feels hot and allow it to cool down thoroughly before using it again.

8.2. Other precautions

Note: Use the initial session to find out how your body responds to vibration therapy. If you experience any discomfort, dizziness or nausea, stop the session immediately.

Note: Do not use the device if any of the contraindications described in chapter 6 apply to you.

Note: Do not use the device for more than 20 minutes per application. Leave it to cool down for at least 15 minutes before restarting.

Note: Make sure to secure the attachments fully onto the novafon device and keep it on during treatment.

Note: If you perform multiple treatment sessions in quick succession with a maximum duration of 20 minutes and only use a single attachment, the temperature of the attachment can reach up to 50 °C. You can avoid this by removing the attachment between sessions (e.g. for cleaning) or by using alternate attachments within the session. Stop the treatment if the contact surface feels uncomfortably warm.

Note: Make sure that no substances, such as liquids, creams or the like, can penetrate into the device or the power plug.

Note: To prevent damage, make sure that the device and charging stand are standing or resting on a stable surface during the charging process.

Note: Do not twist the device in the charging stand. This may cause damage.

Note: Keep the device away from dust, lint and other particles. The penetration of dust, lint or other particles can cause the device to fail. Clean your device regularly with suitable agents to prevent any damage caused by contamination. Please note chapter 9.1.

Note: To prevent premature wear, only store the device in a dry and clean place that is protected against UV light (incl. sunlight).

Note: Use only non-alcoholic cleaning supplies. Alcohol can attack the surfaces of the housing and the attachments.

9. Maintenance and repair

9.1. Cleaning and disinfection

Do not clean the device during charging. Disconnect the device from the power supply when cleaning. If used privately, clean the housing and attachments using a clean cloth moistened with fresh water or a mild soapy solution. For clinical and professional use, clean and disinfect device housings and all attachments with an appropriate surface disinfectant suitable for medical devices before and after each use. As a general rule, only use alcohol-free disinfectants that are suitable for use on sensitive surfaces, such as NOVACLEANER by novafon, SinAlc sensitive by HIMED or RHEOSEPT-WD plus wipes by Rheosol. If necessary, dry off with a soft cloth.

If you are using the charging stand, this can be cleaned with a clean cloth moistened with fresh water or a mild soapy solution if necessary.

To prevent contamination, only place cleaned devices in the charging stand

Ensure that devices are completely dry and with no cleaning agent, disinfectant or oil / cream residues before using the charging stand or charging cable.

9.2. Troubleshooting

If your device automatically switches to a lower setting, or you are unable to increase the intensity, this is probably due to a low battery level. In this case, you should charge your device.

You can also use the FAQs section of our website to find helpful answers to frequently asked questions. This is also available in downloadable format. If other malfunctions should occur in the device, please contact our customer service department.

Never open the device yourself to rectify faults. Any attempt to repair the device yourself will void your warranty.

9.3. Safety check, maintenance and repair

For the novafon device, no regular inspection or maintenance in the sense of a safety check (e.g. according to MPBetreibV) is intended. To ensure safe operation, however, you should inspect the device and attachments for obvious damage such as cracks before use. Attachments must be replaced as soon as they show clear signs of wear or damage (cracks, deformation, changes to the surface). Please note that national laws may nevertheless require you to perform safety checks on the device. Therefore, for professional use, e.g. in practices or hospitals, it is recommended to carry out an electrical safety inspection, e.g. according to DGUV V3 (German Social Accident Insurance), at least every 2 years. This generally does not apply to private individuals. If you experience problems with your device, you may dispose of it or return it to your dealer. Repairs can only be carried out by novafon or qualified dealers. Information on repair options and spare parts is available on our website. Due to cleaning requirements and the necessary IP rating, the battery cannot be replaced by the end user. The battery may only be replaced by novafon or an authorized technician.

Note: Attachments and batteries are wear parts and are therefore subject to limited warranty.

9.4. Disposal

Devices must never be disposed of with common household waste. For environmentally sound disposal, the user is obligated to dispose of the old device at a designated collection point for the disposal of electrical and electronic equipment. The device contains a permanently installed battery and must be disposed of accordingly. Please observe the national

legal regulations when disposing of the device. You can also send your old device to novafon; we will gladly dispose of it for you. This device is subject to EU Directive on Waste Electrical and Electronic Equipment (WEEE). To ensure the conservation of resources, please dispose of the packaging separately as waste paper.

10. Warranty

From the date of delivery, we grant a warranty of 2 years (24 months) for defects resulting from material or production defects. Warranty coverage is subject to the presentation of a proof of purchase. The battery, the attachments is guaranteed for 6 months. Excluded from the warranty: Damage due to improper use, normal wear and tear and consumption (also valid for accessories). The right of rescission (reversal of the purchase), reduction (discount on the purchase price) or subsequent delivery is initially excluded. However, you have a claim to reduction or conversion at your discretion if we do not remedy a defect that is subject to our warranty within a reasonable period of time. Any damage caused by failure to use the product in accordance with its intended use or resulting from improper or unreasonable use or maintenance is not covered by the warranty. Such damage is the sole responsibility of the user. Failure to comply with this instruction manual will void the warranty. Some states may not allow such limitation, so it may not apply to you. This limited warranty gives you specific rights and you may have other rights which may vary from state to state. Get in touch with your dealer or supplier to claim warranty.

11. Customer service

The device may only be repaired by novafon Customer Service. Please do not send in your defective device without being asked, but rather clarify the damage and repair costs with us first. Repair could

be more expensive than a new device. If the repair is no longer financially feasible within the warranty period, novafon offers a replacement unit (again with a 24-month warranty). This manual can also be viewed as a PDF file at www.novafon.com.

12. Simplified EU Declaration of Conformity

The undersigned NOVAFON – Elektromedizinische Geräte GmbH, declares that the product mentioned in this manual complies with the relevant requirements of the following EU regulations:

- Directive 2014/35/EU relating to the making available on the market of electrical equipment designed for use within certain voltage limits
- Directive 2014/30/EU relating to electromagnetic compatibility
- Directive (EU) 2015/863 amending the RoHS Directive
- Regulation (EU) 2023/1542 on batteries and waste batteries
- Regulation (EU) 2025/40 on packaging and packaging waste

The full text of the EU Declaration of Conformity is available at: <https://novafon.com/de/service/downloads#gba>

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