



Release Notes nova fon max+ & max

Exported on 2025-10-01

Versionshinweise - NF Cross-Device-Firmware - 2.11.0 - Mar 26 15:05

Deutsch

- Frequenz und Intensität können jetzt bereits im Welcome-Screen („Hi! FRIEND“) geändert werden
- Verbesserungen bei der Speicherung der Statistiken
- Beim Factory Reset werden gekoppelte Smartphones zuverlässig gelöscht

English

- Frequency and intensity can now be changed already in the welcome screen (“Hi! FRIEND“)
- Improved statistics storage
- During factory reset, paired smartphones are now also deleted – pairing must be repeated afterwards

Versionshinweise - NF Cross-Device-Firmware - 2.10.6 - Oct 24 09:43

Deutsch

- Standard-Intensitätswerte optimiert, damit der Motor noch zuverlässiger startet
- Anpassungen an den Intensitätsstufen für eine bessere Statistik-Auswertung

English

- Optimized default intensity values to ensure the motor starts more reliably
- Adjusted intensity levels for improved statistics tracking

Versionshinweise - NF Cross-Device-Firmware - 2.10.5 - Sep 17 11:24

Deutsch

- Temperaturgrenzen angepasst, um den Antrieb besser zu schützen

English

- Adjusted temperature limits to provide better protection for the drive

Versionshinweise - NF Cross-Device-Firmware - 2.10.4 - Sep 10 14:43

Deutsch

- Initiale Firmware version für NOVAFON max+

English

- Initial Firmware version for NOVAFON max+