

Local vibration therapy  
for animals





1

## PROVEN EFFECT

Pain reduction and improved movement patterns thanks to NOVAFON

The gentle vibrations of the NOVAFON help to reduce pain and alleviate the effects of neurological disorders:



Reduction of tendon pain



Improvement of movement patterns



Relief of muscle pain



Relief of osteoarthritic pain



Relief of joint pain



Promotion of mental and physical relaxation



Improvement of quality of life



Preparation and regeneration of muscles (in sports)



2

## POWERFUL BATTERY

No cables for superior flexibility – at home, at the stables, at the vet or on the go

The new NOVAFON offers complete freedom of movement thanks to its long-lasting battery. Wireless and practical – can be used almost anywhere.

### 3.5 h

Battery life in constant use

### 360°


Freedom of treatment and movement



3

## WIDE RANGE OF TREATMENTS

Intelligent setting options – frequency and intensity can be adjusted depending on the treatment and sensitivity



The new NOVAFON offers 2 frequency settings and 3 intensity settings. The intuitive step-switching enables efficient therapy and the precise setting for the desired treatment mode. From very gentle to strong: The different intensity settings open up a wide range of head-to-toe treatments and can be easily adapted to the individual sensitivity of the animal patient.

## PROVEN TECHNOLOGY

Used successfully for many years

» In every horse that I have treated with the NOVAFON in the test phase, I have noticed a striking relaxation of the muscles treated and a significant warming of the areas treated. «

Animal therapist K. <sup>1,3</sup>



» Dogs with weak muscles in the hindquarters in particular show significant tension in the lower lumbar area. Conversely, permanent pain in this area can also lead to disturbed movement patterns and even to injuries. So the NOVAFON has now found a permanent place in our preparation and warm-up before exercise.«

Veterinary practitioner M. for physiotherapy and naturopathy <sup>2,3</sup>

<sup>1</sup> Homepage, progress report from Hands 4 Horses: [novafon.com/media/pdf/8f/ee/a8/Erfahrungsbericht\\_Pferd\\_Hufe\\_Muskelrelaxation\\_Hands4Horses.pdf](https://novafon.com/media/pdf/8f/ee/a8/Erfahrungsbericht_Pferd_Hufe_Muskelrelaxation_Hands4Horses.pdf). <sup>2</sup> Instagram, therapist report from TierheilPraxis. <sup>3</sup> Individual case studies: There may not be sufficient clinical data to verify the general adequacy of local vibration therapy for the symptoms described. These are quotes from users.

## FOR ANIMAL THERAPISTS AND PET OWNERS

Local vibration therapy for  
horses, cats and dogs

NOVAFON therapeutic treatments, which have been tried and tested on humans for many decades, can also be used on dogs, horses and other animals with orthopaedic or neurological disorders. Targeted local vibrations provide deep stimulation to the muscles, tissues and nerves. This reduces pain and normalises movement patterns triggered by protective postures. Animals perceive this gentle form of therapy as very pleasant and many animal therapists report that, in addition to physiological improvements, animals also show increased psychological relaxation.

## APPLICATIONS:

Pain relief and regeneration in  
various clinical pictures

### HORSES

#### Head, neck and throat muscles

Jaw joint symptoms,  
Neck band symptoms,  
Myofascial pain

#### Shoulder muscles

Myofascial pain

#### Back muscles

Kissing spines,  
Sensitivity in the saddle area,  
Myofascial pain

#### Sacrum joint

Sacroiliac joint symptoms

#### Hindquarters

Knee symptoms,  
Myofascial pain

#### Extremities/legs

Tendon disorders,  
Joint symptoms

### DOGS

#### Head and neck muscles

Myofascial pain

#### Spine and back muscles

Myofascial pain,  
Spondylosis/spondylarthrosis,  
Disc prolapse,  
DLSS – degenerative  
lumbosacral stenosis

#### Extremities

Elbow dysplasia ED,  
Hip dysplasia HD,  
Arthrosis of the joints

### CATS

#### Spine and back muscles

DLSS – degenerative  
lumbosacral stenosis,  
Spondylosis/spondylarthrosis

#### Muscle tension


Myofascial pain

7

## EASY CHANGE OF ATTACHMENTS

The NOVAFON has an intelligent system to mount attachments

All attachments can be easily connected and disconnected by rotating them one quarter turn.



**Disc attachment**  
for extensive use and  
stroking/smoothing



**Ball attachment**  
smaller areas and  
pinpoint treatments

A disc attachment and  
ball attachment are  
always included with  
the device.


8

## FURTHER ATTACHMENTS

Optimise your treatment with  
these attachments\*

### Rehab set

**disc attachment XL** ideal for the  
treatment of larger muscle groups  
**pin attachment** stimulates particularly  
**cone attachment** for selective and  
precise applications



### Magnetic attachment

for very precise and more intense application, ideal for  
treating trigger points and deep-lying areas


### Sensory roller

combines the blood circulation stimulating  
spiky massage ball with vibrations of the  
NOVAFON sound wave device and thus stim-  
ulates the surface and depth sensibility



### Physio set

**fascia ball** suitable for concave surfaces  
**fascia fork** provides punctual pressure, for  
example lateral to the spine  
**fascia wedge** stroking out larger muscle  
strands and deep-lying fascia tissue



\*All attachments are available in white  
and anthracite gray.



9

## PRACTICAL CHARGING STAND

Simplified charging and professional presentation

The charging stand complements the new NOVAFON perfectly: its outstanding features include easier charging thanks to a stable stand and simplified storage of the device. The discreet, timeless design also underlines the device's professional character in the veterinary practice or at home.



10

## INTELLIGENT DETAILS

Various holding options and intuitive one-hand control

The device's ergonomic shape enables users to grip and hold it in a wide variety of ways, thus allowing for efficient, comfortable treatment and targeted application of the vibrations.

All settings can also be controlled intuitively with one hand. This means that the device does not have to be removed mid-treatment, so treatment steps can be seamlessly merged.





Easy to charge  
even without the  
charging stand  
via USB-C















## COMPARISON

The new NOVAFON power –  
The new NOVAFON

Which NOVAFON is best for me?  
You can use our comparison to find out which NOVAFON model is best suited to you and your needs. For the gentler treatment of sensitive areas of the body, we recommend the new NOVAFON. For the treatment of large, compact muscle groups, for example in sport horses, the new NOVAFON power might be the better choice thanks to its wider range of intensities, greater amplitude and more deeply penetrating effect. We are happy to provide you with tailored advice by phone, so that you can find the NOVAFON model best suited to your area of application.

	THE NEW NOVAFON POWER	THE NEW NOVAFON
		
USER GROUPS	Animal therapists and self users for use on sport horses and animals with very compact musculature	Animal therapists and self users for use on cats, dogs and horses (e.g. tendon disorders)
MAXIMUM VIBRATION STRENGTH		
AREAS OF APPLICATION	<p>The entire body – for all clinical pictures</p> <p>Recommended for treatment of compact muscle groups &amp; orthopaedic disorders such as muscle tension, sports injuries and much more.</p>	<p>The entire body – for all clinical pictures</p> <p>Recommended for treatment of sensitive parts of the body &amp; conditions such as tendon disorders, osteoarthritis, scars, CMD and much more.</p>
INTENSITY LEVELS	<p>Intensity 1 </p> <p>Intensity 2 </p> <p>Intensity 3 </p>	<p></p> <p></p> <p></p>
TECHNICAL SPECIFICATIONS		
Battery capacity	3400 mAh	2400 mAh
Battery life	3.5 h	3.5 h
Amplitude	0.3 – 4 mm	0.5 – 3 mm
Penetration effect	6 cm +	6 cm
Treatment levels	100 Hz   50 Hz	100 Hz   50 Hz
Weight of handpiece	245 g	243 g
Dimensions	83 x 187 x 38 mm (W x H x D)	83 x 187 x 38 mm (W x H x D)

## Animal therapists and owners report on their experiences with the NOVAFON

»The NOVAFON relaxes the treated muscle groups, which are blatantly loosened after use. This effect seems to persist over a longer period of time with regular use (two to four times a week over a period of about four weeks) ... Most horses show very clear signs of relaxation when using the NOVAFON: They lower their neck and head, close their eyes and relax the lip muscles.«<sup>4, 6</sup>

Yvonne Katzenberger, Dr. Ruth Katzenberger-Schmelcher,  
PFERGO – 1. Akademie für Pferdeergotherapie

»The NOVAFON is most often used to relieve acute and chronic muscle and joint pain. The gentle vibrations of this sound wave device help reduce pain. But I also use it very successfully for blockages in the musculoskeletal system.«<sup>4, 7</sup>

User N., physiotherapist and naturopath for animals, may 2020,  
Instagram

»Joint problems and muscle tension of a cooperative cat can be treated very well with the NOVAFON and the corresponding magnetic attachment. In any case, Merlin really enjoyed the treatment and commented with a loud purr.«<sup>4, 7</sup>

User A., physiotherapy and osteopathy for animals,  
june 2019, YouTube

<sup>4</sup>Individual case studies: There may not be sufficient clinical data to verify the general adequacy of local vibration therapy for the symptoms described. These are quotes from users. <sup>6</sup>It is a user quote from an experience report that was published on our homepage. <sup>7</sup>It is a user quote that was published on the user's social networks (Facebook or Instagram).

NOVAFON GmbH has been manufacturing electro-medical devices for more than 40 years. With our NOVAFON brand we are an ambassador for local vibration therapy and stand for quality and reliability in this area. Our tried and tested therapy devices are "Made in Germany". As a medical device manufacturer, our quality management is certified according to DIN EN ISO 13485 and our products are subject to strict safety regulations - for the safety and well-being of our customers.



### Certified medical device

NOVAFON is certified as a medical device manufacturer by TÜV Süd. The efficacy of vibration therapy has been investigated and the results have been published in many studies.




### Made in Germany

Our commitment to quality: Highest standards of craftsmanship and material quality. Long-lasting and reliable function. Seamless quality assurance.

Supported by:



on the basis of a decision  
by the German Bundestag

 Products  
[novafon.com/products](https://novafon.com/products)  
Web-Seminars  
[novafon.com/web-seminars](https://novafon.com/web-seminars)

[www.novafon.com](https://www.novafon.com)  
[info@novafon.de](mailto:info@novafon.de)

NOVAFON GmbH  
Daimlerstr. 13  
71384 Weinstadt  
Germany



[novafon\\_vet](#)